
**CompactCath™**

Instruction for use  
Pediatric

LENTIS/MED

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Instructions for use

---

**Product name, model** - CompactCath™, pre lubricated urinary catheter

**Indications for use** - CompactCath™ - intermittent urinary catheter is intended for urological use only. CompactCath is intended to be used by paediatric patients for bladder management, including urine drainage, collection, and urine output measuring. The device is introduced into the urinary bladder via the urethra according to the pictorial instructions on p. 2,3.

**Sterility status** - This product is only sold sterile.

**Storage conditions** - Keep away from sunlight, store at 10°C-25°C, keep dry. After use dispose of the product according to national regulations.

**Contraindications** - No known contraindications.

**Warnings** - This is a single use device. Do not resterilise any portion of this device. Reuse and/or repackaging may create a risk of patient/user infection, compromise the structural integrity and/or essential material and structural characteristics of the device, which may lead to device failure and/or injury, illness or death of the user. Contact a medical specialist in case you notice blood in urine, a reduced quantity or absence of urine, cloudy or thick urine, urine sediment resembling sand, unpleasant smell of urine, frequent urination, high temperature or fever, pain in loins or upper part of abdomen after use.

**Safety measures** - Inspect the catheter before use. Do not use the product if the device or packaging is damaged. Do not use the product after the expiration date.

---

If the application of this product or its use result in a serious harmful event, please report this to the manufacturer and the national supervising body.

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The Lentimed logo is a registered trademark of Lentimed d.o.o.
™ indicates trademarks of Lentimed d.o.o.